
Interface (filters, settings, notifications, sorting, etc.)





This guide offers essential insights into navigating the Platform of therapeutic products interface, highlighting key features such as user settings, notification management, and customizable views.

Version 1.0





Additional settings become visible by selecting the user role at the bottom of the left sidebar.






4 Here you can view and edit the user settings.





Here you can configure notification settings. In-app and email notifications must be actively enabled; otherwise, you will not receive any notifications about new messages, decisions regarding mandatory stock requests, other edits, etc.





In-app notification:

A notification that appears directly within the system – indicated by the bell icon.


Email notification:

Sent to your registered email address – indicated by the envelope icon.

Note: You will only receive notifications (in-app or by email) if you are actively following an issue. An issue can be actively followed by selecting the eye icon on an issue.


By clicking the bell icon at the top left, you can view in-app notifications provided they have been enabled in the notification settings (see previous step).

By clicking on the notification, you will be taken directly to the affected issue.




The **eye icon** indicates that you are following an issue.
If you are the author of an issue, you will automatically follow it.
However, if someone else created the issue, you must **manually activate the eye icon** to receive notifications – or the issue must be **assigned** to you.

8 By clicking on the three dots to the right of the issues overview, you access the view configuration.



In the view configuration, you can select or deselect the columns/properties you want (or do not want) to be displayed.




Issues can be displayed in grouped view – for example, by issue type. When grouping by a specific issue type, **only issues of that type** or **those linked to a issue of that type** will be shown.


Example:

If you group by disruption reports, those will be shown in the main view. Related mandatory stock requests will appear beneath them in an expandable sub-level. If a market overview is linked to the disruption report, it will also be shown – provided you have the required access rights.

However, market requests that are subordinate to the market overview will not be shown.



The grouping function should be used carefully, as not all issues may be visible depending on the filter settings.



You have the option to filter by various dimensions (columns). By selecting "Save as," custom filters can be saved.

**IES NG

| Medicinal products platform | Issues

| Issue | Issue type | Issue type

6/2025

04/06/2025

04/06/2025

21/05/2025

31/05/2025

1 of 3 Pages < >

24/06/2025

18/06/2025

18/06/2025

08/08/2025

31/07/2025

Voriconazol

Buprenorphin-

Rivaroxaban-

HMP-2...

HMP-2...

HMP-2...

HMP-2...

HMP Zulassungsinhaber ... 💠

35t Market request

FBHM.

10 mg, Filmt...

1000 mg, Filmt...

1000 mg, Filmt...

03/06/2025 14:34

03/06/2025 14:10

03/06/2025 11:35


03/06/2025 09:56

03/06/2025 09:25

In Progress

13 Click on "Save" to store the custom filter.

HMP Zulassungsinhaber ... 💠



The newly created filters can be viewed, activated, and modified by selecting the 14 currently activated filter. ***** IES NG ⊚ Test ∨ Sort by Last update ↑ ~ Support Filter New O Test O Test 21/05/2025 27/05/2025 03/06/2025 14:50 03/06/2025 14:46 In Progress 03/06/2025 14:34 Voriconazol 10/06/2025 03/06/2025 14:10 35ug/h, Tra... FBHM, 24/06/2025 HMP-2... Buprenorphin-HMP-2... Rivaroxaban-10 mg, Filmt... FBHM. 04/06/2025 18/06/2025 03/06/2025 11:35 HMP-2... 04/06/2025 18/06/2025 03/06/2025 09:56 35ug/h, Tra... 03/06/2025 09:48 21/05/2025 03/06/2025 09:25 1000 mg, Filmt... 08/08/2025 1000 mg, Filmt... 31/07/2025 03/06/2025 09:25 HMP-2... FBHM, 31/05/2025

1 of 3 Pages < >



Filtered views can be shared via the URL.
The URL contains a **filter ID** that stores your current filter configuration.
By sharing this URL, you can allow other users to view the same specific setup – provided they have the necessary access rights.